Electric Micromotors
Micromotori elettrici
Elektrische Mikromotoren
Micromoteurs électriques
Micromotores eléctricos
Micromotores eléctricos
Mikromotor bezszczotkowy


REF 600.00
REF 602.00

CE
0051

Operating instructions
Istruzioni per l'uso
Betriebsanleitungen
Mode d'emploi
Modo de empleo
Instruções de uso
Instrukcja Obsługi

TKD

ACCESSORIES AND SPARE PARTS / ACCESSORI E RICAMBI / ZUBEHÖR UND ERSATZTEILE / ACCESSOIRES ET RÉCHANGES / ACCESORIOS Y RECAMBOS / ACCESÓRIOS E PEÇAS SOBRESSELENTES / AKCESORIA I CZĘŚCI ZAMIENNE	
REF 001.89.01	Black O-rings / Guarnizioni nere / Schwarze O-Ringe / Joints toriques noirs / Guarniciones negras / O-rings pretos / Czarne O-ringi (10 pcs)
REF 001.98.01	Blue O-rings / Guarnizioni blu / Blauen O-Ringe / Joints toriques bleus / Guarniciones azules / O-rings azules / Niebieskie O-ringi (10 pcs)
REF 098.08.00	LED

TeKne Dental s.r.l.
 Via del Pescinale, 77 - 50041 Calenzano (FI) - Italy
 ☎ +39 055 8825741 ☎ +39 055 8825764
 info@teknedental.com www.teknedental.com

DESCRIPTION

Brushless 3-phase electric micromotors for professional dental use, compliant to IEC 80601-2-60 Standard.
 Micromotors must be employed with any straight or contra-angle handpiece with internal spray and connection compliant to ISO 3964 Standard.
 The DEFINITIVE®LED model is equipped with a LED light source and can be used with any straight or contra-angle handpiece provided with fibre-optics.
Micromotors do not need any lubrication.

This medical device meets the requirements of the European Directives 93/42/EEC (Class IIa) and 2007/47/EC.

- The device:
- must be used by competent and qualified personnel only;
 - is meant exclusively for dental prophylaxis procedures;
 - is not sterilizable;
 - does not produce electromagnetic interference;
 - is not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide;
 - is not liquid-proof (IP 20).

CONTENTS

One DEFINITIVE® or DEFINITIVE®LED micromotor, a set of spare O-rings.

CONNECTION

The DEFINITIVE® (REF 600.00) micromotor must be used together with the specific silicone hose (REF 322.8x).
 The DEFINITIVE®LED (REF 602.00) micromotor must be used together with the specific silicone hose with rotating connector (REF 322.8x.01).
 Insert hose connector into the micromotor rear part (Fig. 1) and screw the metal sleeve (Fig. 2).

Fit the straight or contra-angle handpiece onto the micromotor's coupling and insert it fully.
Never insert the handpiece when the micromotor is rotating.

INSTALLATION

To correctly install or connect the silicone hose to the specific power supply unit (BMC40 electronic board), please refer to the corresponding instructions or application notes.
 The power supply unit must meet the requirements of the CEI EN 60601-1 standard.

On final installation, check dielectric strength and leakage currents according to CEI EN 60601-1 standard.

On installation, micromotor must be cleaned before its first use.

COOLING

Micromotor must be cooled by means of compressed air coming from the dental unit. The minimum recommended air flow is 18 NL/min while the minimum recommended air pressure is 2.0 bar (29 psi).

SPRAY

Maximum spray air pressure: 4.0 bar (58 psi). At 2.0 bar air flow rate is higher than 1.5 l/min.
 Maximum spray water pressure: 4.0 bar (58 psi). At 2.0 bar water flow rate is higher than 50 ml/min.
 The device is equipped with a non-return O-ring for spray water.
 Water must be well-filtered.

CLEANING AND DISINFECTION

After each treatment, clean and disinfect the device immediately. Observing this procedure can easily eliminate any residues of blood, saliva, cooling spray or other.
 The exterior body of the micromotor can be cleaned with a disinfectant containing:
 - glutaraldehyde up to 3% - O-phenylphenol - alcohols based on ethanol
 Disinfectants containing chlorine, acids or solvents are not recommended.
 The device cannot be immersed in any kind of liquid.
The device is not sterilizable.

O.R. REPLACEMENT

Should water or air leakage occur, replacement of O-rings on the coupling is recommended: with the aid of a pin (Fig. 3), pull out the old O-rings and insert the new ones (REF 001.89.01 and REF 001.98.01) on the corresponding grooves.

SERVICE

In case of malfunction of the device, do not try to disassemble micromotor. For any check-up, overhaul or repair work, please contact your supplier, an authorized centre or directly TKD.

GUARANTEE

The device carry a 24-month guarantee against all defects of construction. Guarantee validity is subject to specific conditions.

REFUSE DISPOSAL

This device needs to be recycled. Electrical and electronic equipment may contain dangerous substances which constitute health and environmental hazards. The user must return the equipment to its dealer or establish direct contact with an approved body able to process and derive value from this type of equipment (European Directive 2002/96/EC).

TECHNICAL SPECIFICATIONS

Standard conformity:	ISO 11498: 1997, IEC 80601-2-60: 2012, CEI EN 60601-1: 2007, CEI EN 60601-1-2, European Directives 93/42/EEC and 2007/47/EC.
Classification:	class 2 according to CEI EN 60601-1, type 2 according to ISO 11498 standard
Protection against electrical shock:	type B applied part according to IEC 80601-2-60 standard.
Maximum rotation speed:	40000 rpm
Minimum rotation speed:	depending on the type of used electronic board (1000 rpm with BMC40 board)
Direction of rotation:	clockwise or anti-clockwise
Maximum torque:	depending on the type of used electronic board (3.5 Ncm with BMC40 board)
Maximum peak current:	6 A (10 s)
Maximum noise level:	51 dBA (@45 cm)
Working times:	depending on the type of used electronic board
Cooling air consumption:	26 NL/min @ 3 bar
Sizes:	Ø 21,5 x 77,3 mm
Weight:	90 g
Working temperature:	0 ÷ 40 °C (32 ÷ 104 °F)
Working humidity:	30 ÷ 90 %
Storage temperature:	-10 ÷ 70 °C (14 ÷ 158 °F)
Storage humidity:	10 ÷ 90 %
Light source:	

The LED integrated on the DEFINITIVE®LED micromotor is a source of class 2 according to CEI EN 60825-1:2003 standard. Therefore pay attention not to direct the light beam into eyes.

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