MOON

Operating instructions
Istruzioni per l'uso
Mode d'emploi
Modo de empleo
Instruções de uso
Instrukcja obsługi
Οδηγίες χρήσεως



Fig. 1

REF 520.00	Protection shield / Schermo di protezione / Écran de protection / Pantalla de protección / Escudo de proteção / Oslona zabezpieczająca wzrok / Προστατευτικό καπίρι
REF 530.01	Ø 8 mm Optical fibre / Fibra ottica / Fibre optique / Fibra óptica / Fibra ótica / Światłowód / Οπτική ίνα
REF 530.02	Ø 5 mm Optical fibre / Fibra ottica / Fibre optique / Fibra óptica / Fibra ótica / Światłowód / Οπτική ίνα
<p>ACCESSORIES AND SPARE PARTS / ACCESSORI E RICAMBI / ACCESSOIRES ET RECHANGES / ACCESORIOS Y RECAMBIOS / ACCESÓRIOS E PEÇAS SOBRESSELTANTES / AKCESORIA I CZĘŚCI ZAMIENNE / ΕΣΑΡΤΗΜΑΤΑ ΚΑΙ ΑΝΤΑΓΩΓΙΣΤΙΚΑ</p>	

ENGLISH

DESCRIPTION
Curing lamp for professional dental use.
The lamp can only be used for curing of dental materials suitable for this process (composites, polymer compounds, adhesives, cements) used in various fields of dentistry (traditional, orthodontics, prosthetics).
The light source is provided by a single LED which emits in a narrow band of wavelength, from 430 to 490 nm, with a peak on 460 nm. This band is valid for all materials whose photostarter is made up of camphorquinone (in any case, check whether the material to be cured is sensitive to this wavelength).
To ease operation, lamp can be used with various setups: straight, pistol or in any intermediate position (Fig.2).
The lamp allows 6 curing cycles, different in terms of emitted light intensity and emission times.

Description of the lamp (Fig. 1):

- a) Handpiece
- b) Swivel end part
- c) Optical fibre
- d) Shield to protect eyes
- e) Control keyboard

This medical device meets the requirements of the European Directives 93/42/EEC (Class I) and 2007/47/EC. The device:

- must be used only by trained and qualified personnel;
- has been designed for continuous operation with intermittent load;
- does not produce electromagnetic interference;
- is not suitable for use in presence of a flammable anaesthetic mixture with air, oxygen or nitrous oxide;
- is not liquid-proof (IP20).

CONTENTS

MOON curing lamp, an optical fibre, a protection shield.

CONNECTION

Attach the curing lamp to the connector of the supply silicone hose (REF 330.2x) and screw the corresponding locking nut.
Lamp must be connected to a power supply unit meeting the requirements of CEI EN 60601-1 standard. Particularly, the power supply unit must have a minimum insulation strength of 4000 V between the electric circuits and earth and a minimum insulation strength of 4000 V between the electric circuits and earth.

GENERAL WARNINGS

The LED light is a source of a class 2 according to CEI EN 62471 standard. Such light can damage eyes in case of direct emission without protection.

Never use the curing lamp without the special protection shield on and pay attention not to direct the light beam into eyes.

Patients who suffer from eye pathologies (such as those who have had cataracts removed or who have trouble with the retina) must be adequately protected during lamp use, for example by wearing appropriate protection glasses.
Pay attention to precisely direct the light beam onto the material to be cured. Emitted light can also damage soft tissues (oral mucous membrane, gums, skin).
Do not shake the handpiece excessively and pay attention not to drop it, especially the optical fibre. In case of knock or fall, check handpiece integrity before using the curing lamp. Then turn it on and check its correct operation without using it on a patient.

In case of cracks, or of any other malfunction, do not use lamp on a patient and contact technical service.

OPTICAL FIBRE

To insert the optical fibre onto the handpiece, push it fully down into its housing until a click can be heard.

The optical fibre is made in glass and is therefore fragile: in case of knocks or falls it may break or chip, producing sharp bits which may cause injuries and compromising final quantity of emitted light. If the fibre is damaged, a bright dot appears on the point where it is cracked. In all these cases, the optical fibre must be replaced.

DESCRIPTION OF OPERATING CYCLES

Cycle	Power emission	Energy	Description / LED (Fig.3)
1	1000 mW/cm ² for 20 seconds	20000 mJ/cm ²	Standard cycle. LED 1 turns on.
2	1600 mW/cm ² for 15 seconds	24000 mJ/cm ²	Fast cycle. LED 2 turns on.
3	1800 mW/cm ² for 20 seconds	36000 mJ/cm ²	Strong cycle. LED 3 turns on.
B	500 mW/cm ² for 5 seconds, ramp from 500 to 1000 mW/cm ² for 5 seconds and then 1000 mW/cm ² for 5 seconds.	11250 mJ/cm ²	Bonding. LED 5 + LED 1 turn on. Ramp cycle for a total time of 15 seconds.
R	500 mW/cm ² for 5 seconds, ramp from 500 to 2200 mW/cm ² for 5 seconds and then 2200 mW/cm ² for 5 seconds.	20250 mJ/cm ²	Rapid restoration. LED 5 + LED 2 turn on. Ramp cycle for a total time of 15 seconds.
L	500 mW/cm ² for 5 seconds, ramp from 500 to 1800 mW/cm ² for 5 seconds and then 1800 mW/cm ² for 10 seconds.	26250 mJ/cm ²	Long restoration. LED 5 + LED 3 turn on. Ramp cycle for a total time of 20 seconds.

Data values are valid as output of the @8mm optical fibre.

With the @5mm optical fibre, power and energy values increase by 46 %.

The lamp is equipped with permanent memory so that it powerup the last used cycle is recalled.
The lamp is also equipped with a buzzer which emits a beep when light is turned on, every 5 seconds of operation, and when light is turned off.
Furthermore, the lamp is equipped with a thermal protection.

Sequence of operations

Extract lamp from its holder on the dental unit.

Rotate lamp and/or optical fibre to the best setup for use (straight, pistol or intermediate position).

Select the cycle to be used by means of MODE button.

Put optical fibre into the appropriate position: the output of the optical fibre must be kept as close as possible to the material to be cured, though without touching it. Start the cycle by pressing the START button and wait until light emission turns off automatically.

MODE button (Fig.3)

Press this button to select the cycle to be used. The set cycle is displayed by the turn on of the corresponding LED (or pair of LEDs). Cycle selection is only possible when lamp is not emitting light.

START button

Press this button to start the set cycle. If this button is pressed again during the cycle, light emission will be immediately stopped.

FAULT SIGNALS

In case of malfunction of the curing lamp, the following signals are shown on the key-board:

- LED 5 and LED 1 on: Light beam not being emitted by the lamp. Contact Technical Service.
- LED 5 and LED 2 on: Microcontroller for starting the lamp is malfunctioning. Contact Technical Service.
- LED 5 and LED 3 on: Insufficient power supply. Check supply voltage.
- LED 5 and LED 5 blinking: Handpiece thermal shutdown has tripped. These LEDs will continue to blink until the lamp has cooled down sufficiently (approx. 5 minutes) to be used again. If the problem persists, contact Technical Service.

ENGLISH

CLEANING AND DISINFECTION

Before sterilizing the optical fibre and protection shield, check that there are no residues of cured materials on them; if necessary, remove residues with alcohol or by using a plastic blade.

Curing lamp can be externally disinfected with appropriate products, as for instance germicidal cleansers based on glutaraldehydes, ethanol or propanol.

To disinfect the handpiece use soft disposable paper, avoiding use of corrosive cleansers and avoiding immersing it in any type of liquid.

To disinfect the handpiece, it is advisable to carry out this operation with the optical fibre inserted.

Never use any type of disinfectant on the exposed optical surface of the handpiece when the fibre is removed, as the disinfectant could make this surface become irreversibly cloudy.

STERILIZATION

Both the optical fibre or the protection shield can be sterilized in a steam autoclave at 134 °C/2 bar (273 °F / 29 psi).

Devices are supplied "hot sterile".

The optical fibre can withstand 500 sterilization cycles in autoclave after which it begins to get cloudy and to emit a smaller quantity of light. The protection shield must also be replaced after 500 sterilization cycles.

The handpiece cannot be sterilized in autoclave.

SERVICE

In case of malfunction of the device, and for any overhaul and repair work, please contact your usual supplier or directly to TKD.

REFUSE DISPOSAL

This device needs to be recycled. Electrical and electronic equipment may contain dangerous substances which constitute health and environmental hazards. The user must return the equipment to its dealer or establish direct contact with an approved body able to process and derive value from this type of equipment (European Directive 2002/96/EEC).

GUARANTEE

The device carry a 24-month guarantee against all defects of construction.

TECHNICAL SPECIFICATIONS

- Conformity: ISO 10650-1, CEI EN 60601-1, CEI EN 60601-1-2, Directives 93/42/EEC, 2007/47/EC and 92/31/EEC.
- Degree of protection against direct or indirect contacts: class 1, typeB applied part according to CEI EN 60601-1 standard.
- LED emission wavelength: 430 – 490 nm (peak at 460 nm)
- LED classification: class 2 source, according to CEI EN 62471 standard
- LED lifespan: 5400000 cycles
- Maximum power emission: 2200 mW/cm² ± 10%
- Maximum curable thickness: per single cycle: 3 mm
- Handpiece weight: 120 g
- Supply voltage: 20 ± 36 Vdc or 24 ± 10% Vac, 50/60 Hz; Maximum current consumption: 350 mA
- The power supply unit must meet the requirements of CEI EN 60601-1 standard

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