
SONOSURGERY® STERIL

REF 316.40.08

Rapid coupling with tube
Raccordi rapidi con tubetto
Schnellkupplungen
mit Röhrchen
Raccords rapides avec tube
Enlaces rápidos con tubito

Operating instructions
Istruzioni per l'uso
Bedienungsanleitung
Mode d'emploi
Modo de empleo

Fig. 1

Fig. 2


Fig. 3

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ENGLISH

DESCRIPTION

Surgical coupling for use with the sonic handpiece AIRSON® or SONOSURGERY® or with any other instrument with MULTiflex® connection (Fig. 1). The coupling is intended for driving the handpiece and for the external supply of sterile medical liquids suitable for use in surgery (no saline solution). It is equipped with a structural closing mechanism of the water supply to stop water from the dental treatment unit from flowing through the coupling. Spray air is also stopped.

Maximum air pressure: 6 bar

This medical device meets the requirements of the European Directives 93/42/EEC (Class I) and 2007/47/EC.

CONTENTS

Surgical coupling, tool for installing the coupling, set of spare O-Rings.

CONNECTION TO THE SUPPLY HOSE

Screw the surgical coupling to the supply hose of the dental treatment unit. The connection complies with the standard ISO 9168.

For perfect locking, use the enclosed tool and insert it in the appropriate groove on the coupling (Fig. 2).

To remove the coupling, unscrew the coupling nut of the supply hose and pull off the coupling.

SUPPLY OF STERILE COOLING AGENTS

Make sure to turn off the cooling agent supply of the dental treatment unit. A quantity of at least 50 ml/min of a liquid suitable for medical use (the use of saline solutions is not permitted) must be supplied via external supply systems (components of surgical motors, pressure compressors or syringes). These systems are connected to the lateral connection of the coupling via a supply hose. The constant cooling of the operative site must be guaranteed throughout the entire course of the operation.

LUBRICATION OF THE O-RING

To ensure easy sliding and to improve the impermeability of the coupling, we recommend lubricating the O-rings with small quantities of non-toxic grease in regular intervals.

EXCHANGE OF THE O-RING

We recommend exchanging the O-rings as soon as they leak water or air. The old O-ring can be removed with a pin and the new O-ring can be placed in the corresponding groove (Fig. 3).

CLEANING AND DISINFECTION

The couplings can be cleaned from the outside with non-chlorinated, germicidal agents, provided that these do not contain any Glutaraldehyde, acids or solvents.

Due to the design of the components, an alternative reprocessing of the internal geometries in the bur bath or ultrasonic device is not permitted.

STERILISATION IN THE AUTOCLAVE

For sterilisation, make sure that the packaging is suitable for the coupling which is classed as risk group critical B. The packaging must be large enough to ensure that there is no pressure on the seal.

For sterilisation, the following instructions have to be observed:

Fractionated pre-vacuum (type B)
Sterilisation temperature: 134 °C
Hold time: at least 5 minutes (full cycle)
Drying time: at least 10 minutes

In order to prevent staining and corrosion, the steam must be free of particles (as per DIN EN ISO 13060).

SURGICAL COUPLINGS ARE DELIVERED IN NON-STERILE CONDITION. THEREFORE, THEY HAVE TO BE PREPARED PRIOR TO FIRST USE.

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