**Ultrasonic scaler
Ablatore a ultrasuoni
Ultraschall-Zahnsteinenferner
Détartreur à ultrasons
Detartador ultrasónico**

REF 400.00.01

CE
0051

**Operating instructions
Istruzioni per l'uso
Betriebsanleitungen
Mode d'emploi
Modo de empleo**

Fig. 1


Fig. 2


Fig. 3

Fig. 4

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EMS® is a registered trademark of Electron Medical System SA, Switzerland.

TeKne Dental s.r.l.

Via del Pescinale, 77 - 50041 Calenzano (FI) - Italy
info@teknedental.com www.teknedental.com

DESCRIPTION

Ultrasonic piezo-electric scaler for professional dental use. The scaler can be used for prophylaxis procedures employing the TKD inserts listed below or the EMS® inserts.

This medical device meets the requirements of the European Directives 93/42/EEC (Class IIa) and 2007/47/EC.

The scaler:
- must be used only by trained and qualified personnel;
- is meant for intermittent use and the working times are shown in the technical specifications;
- does not produce electromagnetic interference;
- is not suitable for use with in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide.

CONTENTS

TITANUS®E scaler, a spare O-ring for the handpiece connector, a 24-month guarantee.

CONNECTION OF THE HANDPIECE

After positioning the socket, the handpiece can be fitted to the original EMS® hose by means of a simple snap connection. Do not attempt to screw.

The scaler must be connected to a power supply unit delivering maximum output voltage of 125 V_{rms}. Additionally, the power supply unit must meet the requirements of the CEI EN 60601-1 standard.

COOLING

Operating water pressure for cooling spray can be 0.6 ÷ 2.5 bar. Spray flow rate is higher than 50 ml/min. Cooling spray must be well-filtered.

As inserts heat up when used dry causing heat damage to the teeth, avoid to work dry or, in any case, use inserts without cooling for a few seconds only.

DESCRIPTION OF VIBRATIONS

All ultrasonic vibrations generated by the TITANUS®E transducer move longitudinally along the axis of the insert. Inserts therefore vibrate on a single plane and lateral vibrations, usually common to other instruments, do not exist. For this reason, during procedure, always hold the insert tangentially to the tooth surface.

Do not direct the insert frontally to the enamel. Never direct the tip of the insert vertically to the tooth.

FITTING OF INSERTS

Select the right shape of insert according to the treatment to be performed and screw it manually on to the handpiece (Fig. 1) by means of the dynamometric wrench. Do not fit inserts while the scaler is on.

MAINTENANCE OF THE INSERTS

The shape and weight of each insert are the result of an exact design and manufacturing so as to obtain maximum performance from the ultra-sonic transducer. Therefore any alteration due to fitting, twisting or wearing of the insert can cause anomalous operation: in such case inserts must be systematically replaced. After every procedure, unscrew the insert by means of the dynamometric wrench and sterilise.

USE OF SUPPLIED AND OPTIONAL INSERTS (Fig. 2/3)

INSERT # AX	Use: removal of supragingival tartar deposits on all tooth surfaces. For hard deposits, it is possible to set high power. Power setting (min. ÷ max): 3 ÷ 10
INSERT # BX	Use: removal of supragingival tartar, coatings and stains on all tooth surfaces. Power setting (min. ÷ max): 3 ÷ 10
INSERT # CX	Use: removal of heavy supragingival deposits of tartar. Power setting (min. ÷ max): 3 ÷ 10
INSERT # PX	Use: removal of subgingival tartar deposits on all tooth surfaces, including interproximal and sulcus areas. Power setting (min. ÷ max): 1 ÷ 5
INSERT # PSX	Use: removal of subgingival tartar on all root surfaces, supragingival maintenance therapy. Do not set high power. Power setting (min. ÷ max): 1 ÷ 4
INSERT # DX	Use: removal of crowns, bridges and pivots. Power setting (min. ÷ max): 3 ÷ 8
INSERT # EX	Special 120° ENDO chuck for use with ultrasonic files with 0.8 mm shank diameter. Power setting (min. ÷ max): 2 ÷ 5

The listed power values refer to settings which can be read on the potentiometer dial or LCD display of the power supply unit. Suggested power settings are mean values given as a guide only.

REPLACEMENT OF O-RING

Should water leakage occur, replacement of the O-Ring (REF 001.22.03) on the handpiece connector is recommended. With the aid of a pin (Fig. 4), pull out the old O-ring and insert the new one on the corresponding groove.

STERILIZATION

Handpiece, wrench and inserts can be sterilised in a steam autoclave up to 135 °C / 2 bar (275 °F / 29 psi). Use of a chemiclave sterilizer is prohibited. The handpiece must always be sterilised without the insert on. Use a sterilization bag. The insert instead can be fit in and sterilized together with the dynamometric wrench: this operation avoids any injury or contamination as it is not necessary to touch the insert. Before sterilisation, remove residues and clean the instrument carefully. The device is supplied "not sterile".

CLEANING AND DISINFECTION

Handpiece and wrench can be externally disinfected with a non-chlorinated germicidal cleanser as long as the cleanser does not contain any glutaraldehydes, acid or solvent. Handpiece and wrench cannot be immersed in any kind of disinfectant liquid. For a thorough cleaning of the handpiece, it is possible to unscrew its front cap: this operation allows to disinfect the inner parts, where various residues can build up. Inserts can be cleaned in an ultrasonic tank or with a germicidal cleanser. The screw thread can be properly cleaned with a cotton bud soaked in alcohol.

SERVICE

In case of malfunction of the device and for any overhaul and repair work, please contact your usual supplier or apply directly to TKD.

TECHNICAL SPECIFICATIONS

Standard conformity:	ISO 22374, CEI EN 60601-1, CEI EN 60601-1-2
Operating frequency:	28000 ÷ 31000 Hz.
Classification:	class 1 according to CEI EN 60601-1
Protection against electrical shock:	type B applied part according to CEI EN 60601-1 standard.
Weight without insert:	45 g.
Working temperature:	0 ÷ 40 °C (32 ÷ 104 °F)
Working humidity:	30 ÷ 90 %
Storage temperature:	-20 ÷ 100 °C (-4 ÷ 212 °F)
Storage humidity:	10 ÷ 90 %
Working time with active cooling spray and at minimum power:	maximum 30 min. ON, minimum 10 min. OFF
Working time with active cooling spray and at maximum power:	maximum 8 min. ON, minimum 10 min. OFF
Working time without active cooling spray and at minimum power:	maximum 6 min. ON, minimum 10 min. OFF
Working time without active cooling spray and at maximum power:	maximum 2 min. ON, minimum 10 min. OFF

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